


GEMÜ



**W600 Valve Configurations
in stainless steel**



Leading the world in pharmaceutical and biotechnology industry sterilisation processes

GEMÜ is one of the leading manufacturers of valves, measurement and control systems for sterile applications in the pharmaceutical and biotechnology industries. This position is based on GEMÜ's comprehensive investments in application-oriented research & development, amounting to more than 5% of the company's turnover. The versatile product range is supplemented with a wide range of advisory services provided by industry specialists and application experts.

Customized solutions for your project business

GEMÜ provides the optimal solution from a single source. As a system supplier of isolation, actuator and control technology, we can respond very flexibly to your individual project-specific needs. Our worldwide sales network provides fast reaction times, customer oriented service and a committed project management team.



Table of contents



W600 welding configurations	4	Butt weld spigot bodies	14 - 15
3D and 6D rule	5	Clamp bodies	16
Welding configurations - Selection table	6 - 7	Surface finish	17
GMP / SAP configuration	8	Selection of operators	18 - 19
W600 i-bodies	9	Selection of diaphragms	20
i-bodies - Selection table	10 - 11	Materials and certificates	21
Technical feasibility	12	Specification W600 Valve configurations	22
i-bodies with bypass	13	Worldwide presence	23

W600 welding configurations


The arrangement of two valves welded together to suit the respective application provides maximum functionality in a restricted space. The assembly does without a T piece and thus the dead space between the valves is essentially reduced and two welds are no longer necessary. If more compact designs are required, we recommend using GEMÜ i-bodies and multi-port valve blocks from the GEMÜ M600 series which are machined from a single block. They also have a lower hold-up volume and only a minimum of welds..

Features

- Standard valve body material 1.4435 in investment cast, forged and block material design or 1.4539 in forged design (diaphragm sizes 8 to 50) and block material design
- Various connections selectable
- Various grades of surface finish available
- Operators from the GEMÜ modular system
- Cost effective
- No T piece required
- Valve 2 can be welded on with draining angle



Configuration 2



Configuration 5

3D and 6D rule



Various regulations form the basis for plant designs. Plant operators are normally concerned with the FDA/GMP directives and the ASME/BPE standard. Both regulatory codes define exact geometric reference points for valve configurations. This rule describes the maximum permissible pipe section with a non-turbulent flow in a valve configuration between valve 1 and valve 2. This is either designated as the 3D (3 x dia. ID) rule or the 6D (6 x dia. ID) rule.

3D rule

The longitudinal distance from the main valve inside diameter lower edge to the welded-on sampling valve body sealing weir centre may not exceed 3-times the welded-on sampling valve body inside diameter.

6D rule

The longitudinal distance from the main valve inside diameter centre axis to the welded on sampling valve body sealing weir centre may not exceed 6-times the welded-on sampling valve body inside diameter.




Welding configurations


Selection table

<p>Configuration 1</p>	<p>Configuration 2</p>
<p>Configuration 3</p>	<p>Configuration 4</p>
<p>Configuration 5</p>	<p>Configuration 6</p>
<p>Configuration 7</p>	<p>Configuration 8</p>
<p>Configuration 9</p>	<p>Configuration 10</p>


Configuration 11




Configuration 12



Configuration 13




Configuration 14




Configuration 15

Main body M600, M600 03-01.ER




Configuration 16

Main body M600, M600 03-01.EL




Configuration 17

Main body M600, M600 03-01.ER



Configuration 18

Main body M600, M600 03-01.EL



Notes:

- * Since the max. diameter that can be welded on is limited, we ask that the GEMÜ specification sheet (see page 22) is always used to request the desired combinations
- * The illustrations show recommended installation positions
- * The arrows in the flow charts are examples

S1, S2, S3: Spigots

V1, V2: Valve seat

H: Horizontal orientation

V: Vertical orientation

→ : Flow direction

→ : Draining direction

GMP / SAP configuration




As a rule, the nominal sizes of the two valves differ for GMP and SAP valve configurations. Combinations with the same nominal sizes can, however, also be produced. However, due to the valve geometries and the available space situation (e.g. relating to the actuator dimensions and body), there are also limitations. In these cases, GEMÜ is also able to offer multi-port valve blocks (series M600) manufactured from a single piece as a further customised solution.

SAP valve


The term SAP (Sterile Access Port) valve defines a configuration of two valves welded together, with the 2/2-way valve (1) being arranged horizontally. The valve (2) is welded on vertically in front of or behind the 2/2-way valve (1) sealing weir depending on the application.

GMP valve

The term (Good Manufacturing Practice) valve defines a configuration of two valves welded together, with the 2/2-way valve (1) being arranged vertically. The valve (2) is welded on horizontally in front of or behind the 2/2-way valve (1) sealing weir depending on the application. It is twisted axially to the extent that its sealing weir is turned away from the volumetric flow and that the working medium can flow out unhindered even under depressurised conditions.



SAP valve



GMP valve

W600 i-bodies

The GEMÜ i-body (integrated valve seat) can be seen as an intermediate step to full GEMÜ M-block design machined from a piece of block material. i-bodies are a special construction type of the classical 2/2-way valve bodies. The integrated valve seat of i-bodies is used for example as sampling, steam and condensate valve. The valve bodies have two valve seats and 3 pipe connections. They are manufactured from a forging blank or a piece of block material. The i-body offers a low cost and good alternative for a number of combinations. It already exhibits two essential features of an M-block. It has a greatly reduced dead volume and no internal weld. The drain or supply spigot is only welded on behind the valve seat.

Features

- Reduced weight
- Minimal deadleg
- No weld in the product area
- Compact
- Cost effective
- Available with spigots or elbows
- Draining in vertical installation position possible if adhering to the 3D-rule



i-bodies

Selection table

	IO L	IO R	I1 L	I1 R	I2 L	I2 R
Weld-on parts	None	None	Pipe	Pipe	90° elbow	90° elbow
Flow chart						
Forged bodies						
Block material bodies						

I3L	I3R	I4L	I4R	I5L	I5R
90° elbow	90° elbow	90° elbow	90° elbow	90° elbow	90° elbow

Notes:

- * Alternative installation positions are possible
- * The arrows in the flow charts are examples

S1, S2, S3: Spigots

V1, V2: Valve seat

H: Horizontal orientation



V: Vertical orientation

→ : Flow direction

→ : Draining direction

Technical feasibility

Alternative installation position




Application examples:

- Condensate valve
- Sampling valve

Available seat sizes for material 1.4435:

- Diaphragm size 8/8 block material body
- Diaphragm size 10/8 block material body
- Diaphragm size 25/8 forged body (forged body 1.4539 also possible)
- Diaphragm size 40/8 forged body (forged body 1.4539 also possible)
- Diaphragm size 50/8 forged body (forged body 1.4539 also possible)
- Diaphragm size 80/10 forged body
- Diaphragm size 100/10 forged body



i-bodies with bypass


The bypass valve comprises a forged body with an integrated smaller valve seat. Compared to the i-bodies previously described, this seat is switched in parallel to the main seat. This valve type is particularly suitable for applications which require flow rates that change and that at the same time vary considerably – as is often the case on tank systems and filling machines.

Features


- Variable volumetric flow
- Continuous minimum flow rate
- Precise dosing option
- Large control range thanks to parallel connection of two different seat sizes in one valve
- Works like a static mixer
- Suitable as a sampling and tapping valve
- Different installation positions are possible

Available seat sizes for material 1.4435:

- Diaphragm size 25/8 forged body (forged body 1.4539 also possible)
- Diaphragm size 40/8 forged body (forged body 1.4539 also possible)
- Diaphragm size 50/8 forged body (forged body 1.4539 also possible)
- Diaphragm size 80/10 forged body
- Diaphragm size 100/10 forged body
- Diaphragm size 100/25 forged body
- Diaphragm size 100/40 forged body




Front view



Rear view

Butt weld spigot bodies


The difference between tube specifications (Example DN 15)



Dimensions in mm						DIN		DIN EN 10357 *				DIN 11850		DIN 11866		EN ISO 127 / DIN EN 10357				
MG	DN	NPS	L	LS	H1	Series 0 Code 0	Series B Code 16	Series A Code 17	Series 3 Code 18	Series A Code 1A	Series B Code 1B	Series C Code 60	Series 3 Code 18	Series A Code 1A	Series B Code 1B	Series C Code 60	Series C Code 60			
						ø d	s	ø d	s	ø d	s	ø d	s	ø d	s	ø d	s			
8	4	-	72	20	8.5	6	1.0	-	-	-	-	-	-	-	-	-	-			
	6	-	72	20	8.5	8	1.0	-	-	-	-	-	8	1.0	10.2	1.6	10.2	1.6		
	8	¼"	72	20	8.5	10	1.0	-	-	-	-	-	10	1.0	13.5	1.6	13.5	1.6		
	10	⅜"	72	20	8.5	-	-	12	1.0	13	1.5	14	2.0	13	1.5	-	-	-	-	
	15	½"	72	20	8.5	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
10	10	⅜"	108	25	12.5	-	-	12	1.0	13	1.5	14	2.0	13	1.5	17.2	1.6	17.2	1.6	
	15	½"	108	25	12.5	18	1.5	18	1.0	19	1.5	20	2.0	19	1.5	21.3	1.6	21.3	1.6	
	20	¾"	108	25	12.5	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
25	15	½"	120	25	13.0	19.0	18	1.5	18	1.0	19	1.5	20	2.0	19	1.5	21.3	1.6	21.3	1.6
	20	¾"	120	25	16.0	19.0	22	1.5	22	1.0	23	1.5	24	2.0	23	1.5	26.9	1.6	26.9	1.6
	25	1"	120	25	19.0	19.0	28	1.5	28	1.0	29	1.5	30	2.0	29	1.5	33.7	2.0	33.7	2.0
40	32	1 ¼"	153	25	24.0	26.0	34	1.5	34	1.0	35	1.5	36	2.0	35	1.5	42.4	2.0	42.4	2.0
	40	1 ½"	153	25	26.0	26.0	40	1.5	40	1.0	41	1.5	42	2.0	41	1.5	48.3	2.0	48.3	2.0
50	50	2"	173	30	32.0	32.0	52	1.5	52	1.0	53	1.5	54	2.0	53	1.5	60.3	2.0	60.3	2.0
80	65	2 ½"	216	30	-	62.0	-	-	-	70	2.0	-	-	70	2.0	76.1	2.0	76.1	2.0	
	80	3"	254	30	-	62.0	-	-	-	85	2.0	-	-	85	2.0	88.9	2.3	88.9	2.3	
100	100	4"	305	30	-	76.0	-	-	-	104	2.0	-	-	104	2.0	114.3	2.3	114.3	2.3	

MG = diaphragm size
* replaces DIN 11850

Continued on the next page



Optimum draining angle see brochure „2/2-Way and T Valve Bodies in Stainless Steel“

Dimensions in mm						JIS-G 3447 Code 35	JIS-G 3459 Code 36	SMS 3008 Code 37	BS 4825 Code 55	ASME BPE Code 59	ANSI/ASME B36.19M 10s Code 63	ANSI/ASME B36.19M 40s Code 65								
MG	DN	NPS	L	LS	H1	ød	s	ød	s	ød	s	ød	s	ød	s	ød	s	ød	s	
8	4	-	72	20	8.5	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
	6	-	72	20	8.5	-	-	10.5	1.20	-	-	-	-	-	-	10.3	1.24	10.3	1.73	
	8	¼"	72	20	8.5	-	-	13.8	1.65	-	-	6.35	1.2	6.35	0.89	13.7	1.65	13.7	2.24	
	10	⅜"	72	20	8.5	-	-	-	-	-	-	9.53	1.2	9.53	0.89	-	-	-	-	
	15	½"	72	20	8.5	-	-	-	-	-	-	12.70	1.2	12.70	1.65	-	-	-	-	
10	10	⅜"	108	25	12.5	-	-	17.3	1.65	-	-	9.53	1.2	9.53	0.89	17.1	1.65	17.1	2.31	
	15	½"	108	25	12.5	-	-	21.7	2.10	-	-	12.70	1.2	12.70	1.65	21.3	2.11	21.3	2.77	
	20	¾"	108	25	12.5	-	-	-	-	-	-	19.05	1.2	19.05	1.65	-	-	-	-	
25	15	½"	120	25	13.0	19.0	-	-	21.7	2.10	-	-	-	-	21.3	2.11	21.3	2.77		
	20	¾"	120	25	16.0	19.0	-	-	27.2	2.10	-	-	19.05	1.2	19.05	1.65	26.7	2.11	26.7	2.87
	25	1"	120	25	19.0	19.0	25.4	1.2	34.0	2.80	25.0	1.2	-	-	25.40	1.65	33.4	2.77	33.4	3.38
40	32	1 ¼"	153	25	24.0	26.0	31.8	1.2	42.7	2.80	33.7	1.2	-	-	-	-	42.2	2.77	42.2	3.56
	40	1 ½"	153	25	26.0	26.0	38.1	1.2	48.6	2.80	38.0	1.2	-	-	38.10	1.65	48.3	2.77	48.3	3.68
50	50	2"	173	30	32.0	32.0	50.8	1.5	60.5	2.80	51.0	1.2	-	-	50.80	1.65	60.3	2.77	60.3	3.91
80	65	2 ½"	216	30	-	62.0	63.5	2.0	76.3	3.00	63.5	1.6	-	-	63.50	1.65	73.0	3.05	73.0	5.16
	80	3"	254	30	-	62.0	76.3	2.0	89.1	3.00	76.1	1.6	-	-	76.20	1.65	88.9	3.05	88.9	5.49
100	100	4"	305	30	-	76.0	101.6	2.0	114.3	3.00	101.6	2.0	-	-	101.60	2.11	114.3	3.05	114.3	6.02

MG = diaphragm size


Clamp bodies

All clamp connections are machined according to the spigot dimensions to DIN EN 10357 (replaces DIN 11850), SMS 3008 or ASME BPE. We ask our customers to state which version or standard the connections shall comply with. Depending on the version, clamps are machined from the solid forged body or welded on. Investment cast bodies have welded on clamps as standard.

Welding is carried out by specially qualified and certified welders utilising state-of-the-art welding technology. In principle, special connections requested by customers can be provided on GEMÜ butt weld spigot bodies. Thus it is also possible to have different connections on one body.

Clamp connections for forged 2/2-way bodies	Code
Clamps ASME BPE for pipe ASME BPE, short design	80
Clamps following ASME BPE for pipe EN ISO 1127, length EN 558-1, series 7	82
Clamps ASME BPE for pipe ASME BPE, length EN 558-1, series 7	88
Clamps DIN 32676 for pipe DIN 11850 length EN 558-1, series 7	8A
Clamps SMS 3017 for pipe SMS 3008 length EN 558-1, series 7	8E
Clamps IDF/ISO for pipe JIS-G 3447 length EN 558-1, series 7	8F
Clamps IDF/ISO for pipe JIS-G 3459 length EN 558-1, series 7	8H

Other versions on request



Optimum draining angle see brochure „2/2-Way and T Valve Bodies in Stainless Steel“

Pipe		Code 59 ASME-BPE			Code 60 EN ISO 1127			Code 59 ASME-BPE			Code 16,17,18 DIN 11850			Code 37 SMS 3008			Code 35 JIS-G3447			Code 36 JIS-G3459					
Clamp connection		Code 80			Code 82			Code 88			Code 8A			Code 8E			Code 8F			Code 8H					
DN	NPS	MG	H1	$\phi d1$	$\phi d3$	L	$\phi d1$	$\phi d3$	L	$\phi d1$	$\phi d3$	L	$\phi d1$	$\phi d3$	L	$\phi d1$	$\phi d3$	L	$\phi d1$	$\phi d3$	L	$\phi d1$	$\phi d3$	L	
8	1/4"	8	8	4.57	25	63.5	10.30	25.0	63.5	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10	3/8"		8	7.75	25	63.5	-	-	-	-	-	-	10.00	34	88.9	-	-	-	-	-	-	-	-	-	-
15	1/2"		8	9.40	25	63.5	-	-	-	9.40	25	108	-	-	-	-	-	-	-	-	-	-	-	-	-
10	3/8"	10	12.5	-	-	-	14.00	25.4	108	-	-	-	10.00	34	108	-	-	-	-	-	-	-	14.00	34	108
15	1/2"		12.5	9.40	25	88.9	18.10	50.5	108	9.40	25	108	16.00	34	108	-	-	-	-	-	-	17.50	34	108	
20	3/4"		12.5	15.75	25	101.6	-	-	-	15.75	25	117	-	-	-	-	-	-	-	-	-	-	-	-	-
15	1/2"	25	19	9.40	25	101.6	18.10	50.5	108	9.40	25	108	16.00	34	108	-	-	-	-	-	-	17.50	34	108	
20	3/4"		19	15.75	25	101.6	23.70	50.5	117	15.75	25	117	20.00	34	117	-	-	-	-	-	-	-	-	-	
25	1"		19	22.10	50.5	114.3	29.70	50.5	127	22.10	50.5	127	26.00	50.5	127	22.60	50.5	127	23.00	50.5	127	-	-	-	
32	1 1/4"	40	26	-	-	-	38.40	64	146	-	-	-	32.00	50.5	146	31.30	50.5	146	29.40	50.5	146	-	-	-	
40	1 1/2"		26	34.80	50.5	139.7	44.30	64	159	34.80	50.5	159	38.00	50.5	159	35.60	50.5	159	35.70	50.5	159	-	-	-	
50	2"	50	32	47.50	64	158.75	56.30	77.5	190	47.50	64	190	50.00	64	190	48.60	64	190	47.80	64	190	-	-	-	
65	2 1/2"	80	62	60.20	77.5	193.68	72.10	91	216	60.20	77.5	216	66.00	91	216	60.30	77.5	216	59.50	77.5	216	-	-	-	
80	3"		62	72.90	91	222.25	84.30	106	254	72.90	91	254	81.00	106	254	72.90	91	254	72.30	91	254	-	-	-	
100	4"	100	76	97.38	119	292.1	109.70	144.5	305	97.38	119	305	100.00	119	305	97.60	119	305	97.60	119	305	-	-	-	

Dimensions in mm
MG = diaphragm size

Surface finish









Modern, ergonomically shaped workstations and trained polishing staff give us the ability to provide high quality surface finishes. Depending on the required application, surface finishes from Ra 0.8 µm down to 0.25 µm can be achieved by polishing, electro polishing or a special process, we call "elysieren".

Mechanical hand polishing is carried out at our works to ensure our high quality standard. In principle, special connections requested by customers can be provided on GEMÜ butt weld spigot bodies and it is also possible to have different connections on one body.

Valve body surface finish, internal contour			
	Forged body - Codes 40, 42, F4 Block material - Codes 41, 43	Investment casting Codes 32, 34	Code
Ra ≤ 0.8 µm for media wetted surfaces, mechanically polished internal	X	X	1502
Ra ≤ 0.8 µm for media wetted surfaces, electropolished internal/external	X	-	1503
Ra ≤ 0.6 µm for media wetted surfaces, mechanically polished internal	X ¹	X ¹	1507
Ra ≤ 0.6 µm for media wetted surfaces, electropolished internal/external	X ¹	-	1508
Ra ≤ 0.25 µm for media wetted surfaces, electropolished internal/external	X ¹	-	1516
Ra ≤ 0.25 µm for media wetted surfaces, mechanically polished internal	X ¹	-	1527
Ra ≤ 0.4 µm for media wetted surfaces, mechanically polished internal	X ¹	-	1536
Ra ≤ 0.4 µm for media wetted surfaces, electropolished internal/external	X ¹	-	1537
Ra ≤ 0.51 µm (20 µinch) for media wetted surfaces, mechanically polished internal	X ¹	-	1927
Ra ≤ 0.51 µm (20 µinch) for media wetted surfaces, electropolished internal/external	X ¹	-	1928
Ra ≤ 0.38 µm (15 µinch) for media wetted surfaces, electropolished internal/external	X ¹	-	1929

Ra in accordance with DIN 4768; at defined reference points. ¹For pipe inside diameter < 6 mm, surface in spigots Ra ≤ 0.8 µm.

Selection of operators

	Manually operated						Motorized	
								
Type	9601**	9602**	9612**	9673	9653**	9654**	9618	9698
Material	Stainless steel, plastic handwheel, with optical position indicator and seal adjuster	Stainless steel, with optical position indicator and seal adjuster	Stainless steel, plastic handwheel, with optical position indicator and seal adjuster	Stainless steel, plastic handwheel, with optical position indicator and seal adjuster	Stainless steel, plastic handwheel, with optical position indicator, stroke limiter/ seal adjuster, lockable, optional: electrical position indicator	Stainless steel, with optical position indicator, stroke limiter/ seal adjuster, lockable, optional: electrical position indicator	Plastic, with or without stainless steel distance piece, optical position indicator	Plastic, with or without stainless steel distance piece, optical position indicator and manual override
Autoclavable	●	●	●	●	●	●	-	-
Operating temperature*	-10 to 150 °C	-10 to 150 °C	-10 to 150 °C	-10 to 150 °C	-10 to 150 °C	-10 to 150 °C	0 to 130 °C (without distance piece 15 to 60 °C)	0 to 150 °C
Operating pressure*	0 to 10 bar	0 to 10 bar	0 to 10 bar	0 to 10 bar	0 to 10 bar	0 to 10 bar	0 to 6 bar	0 to 6 bar
DN	4 to 15	4 to 15	10 to 20	15 to 50	10 to 100	4 to 100	4 to 15	15 to 50
Supply voltage	-	-	-	-	-	-	24 VAC, 120 VAC, 230 VAC, 50/60Hz	24 VAC, 120 VAC, 230 VAC, 50/60Hz
Diaphragm size 8	●	●	-	-	-	●	●	-
Diaphragm size 10	-	-	●	-	●	●	●	-
Diaphragm size 25	-	-	-	●	●	●	-	●
Diaphragm size 40	-	-	-	●	●	●	-	●
Diaphragm size 50	-	-	-	●	●	●	-	●
Diaphragm size 80	-	-	-	-	●	●	-	-
Diaphragm size 100	-	-	-	-	●	●	-	-

* dependent on diaphragm material, see technical datasheet

** also suitable for i-body for valve seat 2





Pneumatically operated



Type	9605	9625	9687	9650**	9658/9688	9660
Material	Plastic, with stainless steel distance piece, optical position indicator	Plastic, with stainless steel distance piece, optical position indicator	Plastic, with stainless steel distance piece	Stainless steel, with optical position indicator, optionally autoclavable	Two stage actuator, stainless steel	Filling valve, stainless steel with optical position indicator
Autoclavable	-	-	-	(DN 4 to 25)	-	-
Operating temperature*	-10 to 150 °C	-10 to 150 °C	-10 to 150 °C	-10 to 150 °C	-10 to 150 °C	-10 to 150 °C
Operating pressure*	0 to 8 bar	0 to 6 bar	0 to 10 bar	0 to 10 bar	0 to 10 bar	0 to 5 bar
DN	4 to 15	10 to 20	10 to 100	4 to 100	10 to 50	4 to 25
Supply voltage	-	-	-	-	-	-
Diaphragm size 8	●	-	-	●	●	●
Diaphragm size 10	-	●	●	●	●	●
Diaphragm size 25	-	-	●	●	●	●
Diaphragm size 40	-	-	●	●	●	-
Diaphragm size 50	-	-	●	●	●	-
Diaphragm size 80	-	-	●	●	-	-
Diaphragm size 100	-	-	●	●	-	-

* dependent on diaphragm material, see technical datasheet

** also suitable for i-body for valve seat 2



Selection of diaphragms

Diaphragm	Material/Design	Dia-phragm size	Temperature range [°C]			Code	Certificates and approvals				
			Min.	Max.	Sterilisation ¹		FDA compli-ant	USP Class VI	EHEDG	TA Luft (German Clean Air Act)	O ₂ BAM
EPDM	Ethylene-propylene-diene rubber	8 - 100	-10	100	max. 150 °C ² max. 60 min. per cycle	13/3A	☒	☒	☒	☒	☒
EPDM	Ethylene-propylene-diene rubber	8 - 100	-10	100	max. 150 °C ² max. 180 min. per cycle	17	☒	☒	☒	☒	
PTFE/EPDM	Fully laminated PTFE diaphragm with EPDM back	8, 10, 100	-10	100	max. 150 °C ² , no time limit per cycle	52/5A	☒	☒	☒	☒	☒
	Convex two-piece PTFE diaphragm with loose EPDM back	25, 40, 50, 80	-10	100	max. 150 °C ² , no time limit per cycle	5E	☒	☒	☒	☒	☒

¹ The sterilisation temperature is valid for steam (saturated steam) or superheated water.


² If the sterilisation temperatures listed above are applied to the EPDM diaphragms for longer periods of time, the service life of the diaphragms will be reduced. In these cases, maintenance cycles must be adapted accordingly.

This also applies to PTFE diaphragms exposed to high temperature fluctuations.


PTFE diaphragms can also be used as moisture barriers; however, this will reduce their service life.

The maintenance cycles must be adapted accordingly.


GEMÜ 555 and 505 globe valves are particularly suitable for use in the area of steam generation and distribution. The following valve arrangement for interfaces between steam pipes and process pipes has proven itself over time: A globe valve for shutting off steam pipes and a diaphragm valve as an interface to the process pipes.



Das Original GEMÜ Dichtsystem



GEMÜ seal system




Conventional seal systems

Materials and certificates

The table below provides an overview of the possible certificates which are generally available. The type of certificate and its content must be specified exactly before ordering to be able to provide the required documents. Later requests of certificates may not be possible or possible only under certain conditions.

Our specialists are happy to answer any questions you might have.

Type	Designation of the test certificate in accordance with EN 10204	Content of the certificate	Confirmation of the certificate by
2.1	Certificate of compliance with the order	Confirmation of compliance with the order	the manufacturer
2.2	Test report	Confirmation of compliance with the order with specification of results of non-specific testing	the manufacturer
3.1	Inspection certificate 3.1	Confirmation of compliance with the order with specification of results of specific testing	the manufacturer acceptance officer independent of the production division
3.2	Inspection certificate 3.2	Confirmation of compliance with the order with specification of results of specific testing	the manufacturer acceptance officer independent of the production division and the acceptance officer commissioned by the purchaser or the acceptance officer named in the official regulations



Specification W600 Valve configurations

Please complete this form and return it to your nearest GEMÜ office or to the address listed below!

Configuration no.:

Operating pressure

bar

Quantity

Medium temperature

°C

Valve 1

	DN	Code	ød(a) [mm]	s [mm]
Spigot S1				
Spigot S2				

only bored (make entries for dimensions for S3)

Valve 2

	DN	Code	ød(a) [mm]	s [mm]
Spigot S3				

no deadleg requirement

3D rule
See sketch below

6D rule
See sketch below

Operator type

Operator type

Diaphragm size

Diaphragm size

Control function

Control function

Accessories

Accessories

Remark

Remark



Body material * Valve 1	1.4435	<input type="radio"/>
	1.4435 BN 2 (Δ Fe < 0,5%)	<input type="radio"/>
	1.4539	<input type="radio"/>
	Other	<input type="radio"/>
	* Forged body as standard	

Body material * Valve 2	1.4435	<input type="radio"/>
	1.4435 BN 2 (Δ Fe < 0,5%)	<input type="radio"/>
	1.4539	<input type="radio"/>
	Other	<input type="radio"/>
	* Forged body as standard	

Diaphragm material	EPDM	<input type="radio"/>	Code	_____
	PTFE	<input type="radio"/>	Code	_____
	Other	<input type="radio"/>		_____

Diaphragm material	EPDM	<input type="radio"/>	Code	_____
	PTFE	<input type="radio"/>	Code	_____
	Other	<input type="radio"/>		_____

Surface finish internal contour valve 1 and 2	1502	(Ra) ≤ 0,8 µm	<input type="radio"/>
	1503	(Ra) ≤ 0,8 µm electropolished	<input type="radio"/>
	1507	(Ra) ≤ 0,6 µm	<input type="radio"/>
	1508	(Ra) ≤ 0,6 µm electropolished	<input type="radio"/>
	1536	(Ra) ≤ 0,4 µm	<input type="radio"/>
	1537	(Ra) ≤ 0,4 µm electropolished	<input type="radio"/>
	1527	(Ra) ≤ 0,25 µm	<input type="radio"/>
1516	(Ra) ≤ 0,25 µm electropolished	<input type="radio"/>	
Other	_____	<input type="radio"/>	



The technical details of each enquiry must be checked by GEMÜ.

Contact (GEMÜ):	_____
Customer:	_____
Department:	_____
Address:	_____
Phone:	_____
E-mail:	_____

For internal notes

Worldwide presence

AUSTRALIA

GEMÜ Australia Pty. Ltd
Unit 4 - 8/10 Yandina Road
West Gosford, NSW 2250
Phone: +61-2-43 23 44 93
Fax: +61-2-43 23 44 96
mail@gemu.com.au

AUSTRIA

GEMÜ GmbH
Europaring F15 401
2345 Brunn am Gebirge
Phone: +43 2236 30 43 45-0
Fax: +43 2236 30 43 45-31
info@gemue.at

BELGIUM

GEMÜ Valves bvba/sprl
Koning Albert 1 laan, 64
1780 Wemmel
Phone: +32 2 702 09 00
Fax: +32 2 705 55 03
info@gemue.be

BRAZIL / SOUTH AMERICA

GEMÜ Indústria de Produtos Plásticos e Metalúrgicos Ltda.
Rue Marechal Hermes, 1141
83.065-000 São José dos Pinhais Paraná
Phone: +55-41-33 82 24 25
Fax: +55-41-33 82 35 31
gemu@gemue.com.br

CANADA

GEMÜ Valves Canada Inc.
2572 Daniel-Johnson Boulevard
Laval, Quebec
H7T 2R8
Phone: +1-450-902-2690
Fax: +1-404-3 44 4003
info@gemu.com

CHINA

GEMÜ Valves (China) Co., Ltd
No.518, North Hengshahe Road
Minhang District, 201108
Shanghai
Phone: +86-21-2409 9878
info@gemue.com.cn

DENMARK

GEMÜ ApS
Industriparken 16-18
2750 Ballerup
Phone: +45 70 222 516
info@gemue.dk

FRANCE

GEMÜ S.A.S
1 Rue Jean Bugatti
CS 99308 Duppigheim
67129 Molsheim Cedex
Phone: +33-3 88 48 21 00
info@gemu.fr

INTERCARAT

1 Rue Jean Bugatti
CS 99308 Duppigheim
67129 Molsheim Cedex
Phone: +33-3 88 48 21 20
sales@intercarat.com

GERMANY

GEMÜ Gebr. Müller GmbH & Co. KG
Fritz-Müller-Straße 6 - 8
74653 Ingelfingen-Criesbach
Postfach 30
74665 Ingelfingen-Criesbach
Phone: +49 (0)7940-12 30
Fax: +49 (0)7940-12 31 92
(Domestic)
Fax: +49 (0)7940-12 32 24 (Export)
info@gemue.de

Inevvo solutions GmbH & Co. KG
Fritz-Müller-Platz 1
74676 Niedernhall-Waldzimmern
Phone: +49 (0)7940-12 38 681
info@inevvo-solutions.com

GREAT BRITAIN

GEMÜ Valves Ltd.
10 Olympic Way
Birchwood, Warrington
WA2 0YL
Phone: +44-19 25-82 40 44
Fax: +44-19 25-82 80 02
info@gemu.co.uk

HONG KONG

GEMÜ (Hong Kong) Co., Ltd.
Room 2015, Tower B,
Regent Centre,
70 TA Chuen Ping Street
Kwai Chung, N.T., Hong Kong
P.R. China
Phone: +852 6873 8280
Fax: +852 6873 8280
info@gemue.com.cn

INDIA

GEMÜ Branch Office India
Room Number 135,
1st Floor, 101-104, B-Wing,
GCP Business Centre
Opp. Memnagar Fire Station
Vijay Cross Road
Ahmedabad - 380 014
Phone: +91-79-6134 4423
Fax: +91-79-25450439
sales@gemu.in

INDONESIA

GEMÜ Valves Pte Ltd
(Indonesia Representative Office)
Rukan Mangga Dua Square
Block F17, 2nd Floor
Jl. Gunung Sahari Raya No. 1
Jakarta Utara 14420
Indonesia
Phone: +62 (21) - 6231 0035
Fax +62 (21) - 2907 4643
info@gemu.co.id

IRELAND

GEMÜ Ireland Ltd
15 Eastgate Drive
Eastgate Business Park
Little Island
Co. Cork
Phone: +353 (0)21 4232023
Fax: +353 (0)21 4232024
info@gemu.ie

ITALY

GEMÜ S.r.l.
Via Giovanni Gentile, 3
20157 Milano
Phone: +39-02-40044080
Fax: +39-02-40044081
info@gemue.it

JAPAN

GEMÜ Japan Co., Ltd.
2-5-6, Aoi, Higashi-ku,
Nagoya, Aichi, 461-0004
Phone: +81-52-936-2311
Fax: +81-52-936-2312
info@gemu.jp

MALAYSIA

GEMÜ Valves PTE. LTD.
D-8-09, Block D, No. 2A
Jalan PJU 1A/7A
Oasis Square, Oasis Damansara
47301 Petaling Jaya
Selangor Darul Ehsan
Phone: +(603)- 7832 7640
Fax: +(603)- 7832 7649
info@gemu.com.sg

MEXICO

GEMÜ Valvulas S.A. de C.V.
German Centre,
Av. Santa Fe No. 170 - OF. 5-1-05
Col. Lomas de Santa Fe,
Del. Alvaro Obregon
01210 Mexico, D.F.
Phone: +52 55 7090 4161
+52 55 7090 4179

RUSSIA

OOO „GEMÜ GmbH“
Uliza Shipilovskaya, 28A
115563, Moskau
Phone: +7(495) 662-58-35
Fax: +7(495) 662-58-35
info@gemue.ru

SINGAPORE

GEMÜ Valves PTE. LTD.
25 International Business Park
German Centre #03-73/75
Singapore 609916
Phone: +65-65 62 76 40
Fax: +65-65 62 76 49
info@gemu.com.sg

SOUTH AFRICA

GEMÜ Valves Africa Pty. Ltd
Cnr Olympic Duel Avenue
And Angus Crescent,
Northlands Business Park
(Stand 379),
New Market Road
Randburg
Phone: +27 11 462 7795
Fax: +27 11 462 4226
info@gemue.co.za

SWEDEN

GEMÜ Armatör AB
Heljesvägen 8
437 36 Lindome
Phone: +46-31-996500
order@gemu.se

SWITZERLAND

GEMÜ GmbH
Seetalstr. 210
6032 Emmen
Phone: +41-41-7 99 05 05
Fax: +41-41-7 99 05 85
info@gemue.ch

GEMÜ Vertriebs AG
Lettenstrasse 3
6343 Rotkreuz
Phone: +41-41-7 99 05 55
Fax: +41-41-7 99 05 85
vertriebsag@gemue.ch

TAIWAN

GEMÜ Taiwan Ltd.
9F.-5, No.8, Ziqiang S. Rd.
Zhubei City
Hsinchu County 302,
Taiwan (R.O.C.)
Phone: +886-3-550-7265
Fax: +886-3-550-7201
office@gemue.tw

UNITED STATES

GEMÜ Valves Inc.
3800 Camp Creek Parkway
Suite 120, Building 2600
Atlanta, Georgia 30331
Phone: +1-678-5 53 34 00
Fax: +1-404-3 44 93 50
info@gemu.com

In addition to these subsidiaries, GEMÜ has a global partner network.

Contact details:
www.gemu-group.com/en_GB/kontakte/



 GEMÜ manufacturing site

 GEMÜ subsidiary

